

# Determination of flame-retardant materials in plastics using a combination of analytical techniques

By Anthony Fong, Dr. Jesse Allen, and Dr. Yanika Schneider, EAG Laboratories

## ABSTRACT

Flame retardant compounds serve an important purpose in society and are particularly critical in plastics, which are often more flammable than other materials. To evaluate the efficacy of flame retardants in commercial products, it is important to know both the concentration and composition. However, due to the variability of flame retardants, the appropriate analytical method is not always obvious. In this publication, we analyze an unknown plastic box advertised to have flame retardant properties. We use a series of analytical techniques and evaluate their compatibility with one another.

## INTRODUCTION

Flame retardant materials constitute an important class of materials for protecting life and property. These compounds can help slow or stop the spread of fire, and as such are typically found in a variety of consumer goods. The exact composition and concentration of these compounds is essential for determining their efficacy, especially in cases of litigation.<sup>1</sup>

Many different retardant materials are commercially available, including brominated, phosphorus, nitrogen, chlorinated, and inorganic compounds such as hydrated aluminum, magnesium oxide, and antimony trioxide.<sup>2-4</sup>

Due to the variety of flame retardant materials, several analytical tools must be employed to determine the composition. In this publication, a commercially available plastic box was analyzed using a combination of fourier transform infrared spectroscopy (FTIR), x-ray fluorescence spectroscopy (XRF) and gas chromatography mass spectroscopy (GCMS) to illustrate the power of complimentary techniques to determine the composition and concentration of the flame retardant component.

## EXPERIMENTAL

A plastic box shown in Figure 1 was purchased from an online



**Figure 1. Photo of flame retardant-filled plastic box.**

retailer. The box was cut into three pieces using a band saw. Experiments were performed using FTIR, XRF and GCMS techniques as described below.

**FTIR:** A small portion of the sample was transferred to an infrared transmitting substrate and examined using a Thermo-Nicolet 6700 Fourier Transform Infrared (FTIR) spectrometer equipped with a Continuum microscope in transmission mode. The analytical spot size was approximately 100 microns x 100 microns. OMNIC 8.0 software was used to perform data analysis.


**XRF:** X-ray Fluorescence (XRF) was performed using a Rigaku Primus II WDXRF with a rhodium X-ray source x-ray tube, vacuum atmosphere and an analysis area of 30 mm diameter. This analysis utilized a wavelength dispersive spectrometer (WDXRF) that is capable of detecting elements from atomic number (Z) of 4 (beryllium) through atomic number 92 (uranium) at concentrations from the low parts per million (ppm) range up to 100% by weight.

**GCMS:** An extract was prepared by dissolving a small piece of plastic in dichloromethane. The extract was analyzed on an Agilent 6890A Gas Chromatograph/ Agilent 5973 Mass Spectrometer with a mass selective detector (MSD) using a 30 m X 0.25 mm DB-5MS (J&W Scientific) column under a flow rate of 1 mL/min. Once the flame retardant material was determined, a 1-point calibration curve was prepared using a commercial standard with a known concentration.

# Determination of flame-retardant materials in plastics using a combination of analytical techniques

## RESULTS AND DISCUSSION

A commercially available plastic box shown in Figure 1 was examined by FTIR in transmission mode. The resulting spectrum is shown in Figure 2 (black box). The spectrum is dominated by peaks at 3062, 3027, 2922, 2848, 2238, 1602, 1475, 1452, 759, and 701  $\text{cm}^{-1}$ , which are due to acrylonitrile butadiene styrene (ABS) as determined by comparison with library references. In addition to ABS, several other peaks were observed at 1551, 1475, 1393, 1318, 278, 1243, 1159, and 737  $\text{cm}^{-1}$ , which are due to tetrabromobisphenol A (TBBPA). TBBPA is a commonly used fire retardant in a variety of applications.<sup>3,4</sup> It is also an initiator polymerization of ABS. Furthermore, a trace ester is observed at 1728  $\text{cm}^{-1}$ , which may be an additive in the plastic. Altogether, these results demonstrate that TBBPA is the main flame retardant component in the plastic and that the base polymer is ABS.




**Figure 2: FTIR spectrum of flame retardant-filled box demonstrating presence of ABS and TBBPA**

XRF was performed directly on a small piece of plastic and the results are summarized in Table 1. A representative XRF spectrum is shown in Figure 3. The main elements detected include C, Br, Sb, N and O. The presence of C and N confirms the FTIR finding of ABS base polymer, while the high bromine level supports the presence of TBBPA as a flame retardant component. The ~10% antimony content also suggests that antimony trioxide may be a second flame retardant component in the plastic, which is also supported by the ratio of Sb to O. Antimony trioxide is an inorganic flame retardant commonly used in combination with brominated compounds.<sup>5</sup>

**Table 1: Elemental composition determined by XRF<sup>a</sup>**

Element	Conc.	Element	Conc.
C	55.3	Mg	0.028
Br	25.0	Na	0.008
Sb	9.73	K	0.008
N	5.59	Ca	0.006
O	3.90	Fe	0.006 <sup>b</sup>
S	0.120	Cl	0.004
F	0.100	Ni	0.002 <sup>b</sup>
Si	0.090	Zn	0.001
P	0.034	Cu	0.001

- a The results are normalized to 100% of the measured and detected elements
- b At levels observed, the concentrations may have a considerable contribution from instrumental background



**Figure 3: XRF spectrum of flame retardant-filled box showing high levels of Br and Sb**

GCMS analysis was performed on a dichloromethane polymer extract and the major volatiles are summarized in Table 2. The presence of TBBPA is again confirmed by the mass spectrum shown in Figure 4. Other compounds are observed in the plastic, including antioxidant(s) and ABS-derivatives. Note that tribromobisphenol A is a de-bromination product of TBBPA.<sup>6</sup>

The concentration of TBBPA in the extract solution was determined using a 1-pt calibration with a known concentration of TBBPA. The amount of TBBPA in the extract was found to be 6.83 mg/mL, which equals approximately 48 wt%. The amount of Br in the extract is then calculated to be approximately 28 wt%. This value is very close to the 25% wt determined by XRF, suggesting that the two techniques provide complimentary results.

# Determination of flame-retardant materials in plastics using a combination of analytical techniques

Table 2: Volatiles detected in polymer extract by GCMS

Avg. RT	Normalized Area/mg	Compound Identification
11.27	0.32	Possibly an antioxidant
11.64	0.17	Antioxidant group
22.59-22.76	0.57	2-[1-(4-Cyano-1,2,3,4-tetrahydronaphthyl)] propanenitrile
29.29	0.45	Tribromobisphenol A
31.38	47.9	Tetrabromobisphenol A




Figure 4: Mass spectrum of polymer extract confirming presence of TBBPA

## CONCLUSIONS

A plastic box containing an unknown flame retardant was examined using a suite of analytical techniques including FTIR, XRF and GCMS. The flame retardant components were identified as tetrabromobisphenol A (TBBPA) and possibly antimony trioxide. The amount of TBBPA was determined to be ~25 wt% by XRF and ~28% by GCMS. The proximity of the two values demonstrates the compatibility of the methods. Altogether, the results suggest that the combination of techniques produces complimentary and reliable information about the flame retardant composition.

## REFERENCES

1. E. D. Weil and S. Levchik, *Journal of Fire Sciences* **22**, 25 (2004).
2. S. Zhang and A.R. Horrocks, *Progress in Polymer Science*, **28**, 1517 (2003).
3. L. Chen and Y.-Z. Wang, *Polym. Adv. Technol.*, **21** 1–26 (2010).
4. E. D. Weil, *Journal of Fire Sciences* **29**, 259 (2011).
5. U. Sellström, B. Jansson, *Chemosphere*, **31**, 3085 (1995).
6. J. Eriksson, S. Rahm, N. Green, Å. Bergman and E. Jakobsson. *Chemosphere*, **54**, 117 (2004).